

KytoMatrix Radial Compression Band Clinical Evidence Summary

Introduction

The KytoMatrix Radial Compression Band is a safe, effective, and rapid patent hemostasis-achieving device used to control bleeding following transradial catheterizations. Through a unique combination of inflatable radial band technology with a built in Chitosan impregnated patch, accelerated blood clotting¹⁻³ and infection reducing capabilities are achieved.⁴ The quicker bleeding control is independent of the blood-clotting cascade and is complementary to normal clotting, allowing for effective and safe bleeding control in individuals with bleeding disorders, patients on anticoagulant medication, and children.⁵⁻⁷

The more efficient control of bleeding decreases the length of stay and time to discharge, owing to the possibility of an increase in procedural cases being scheduled.⁸ A direct correlation exists between a reduction in patient complications and the use of radial access combined with radial band technology; this is facilitated by patent hemostasis achieved with shorter compression times and low compression pressures.^{1,9-13}

The approximate use of oral anticoagulants in treatment visits is more than 8 million per year.⁴⁶

Accelerate the recommended times

Continuous firm pressure should be held at the site to obliterate the distal pulse for **3-5 minutes**. After the initial 3-5 minutes pressure must be reduced enough for the healthcare provider to obtain a distal pulse.⁴⁴

- General recommendation is **4 minutes** per French
- Hold pressure for at least **15-20 minutes** for a diagnostic procedure
- Hold pressure for at least **30 minutes** for an interventional procedure
- Hold pressure for at least **45 minutes** for the removal of balloon pump or larger bore sheaths

Suggested cardiac catheterization protocols using the KytoMatrix™ Dressing (Chitosan impregnated patch).⁴⁵

Compression times diagnostics (5F-8F)	5-10 minutes
Compression times interventional (5F-8F)	10-15 minutes

Fewer Complications and Better Patient Experience

The radial band is designed with key features that enhance patient outcomes with respect to the incidence of radial artery occlusion (RAO) and hematoma rates as compared to manual compression.^{14,15} This is evident through **a statistically significant reduction in RAO and hematoma occurrence with use of band technology.**^{10,14,16}

The combination of a band with different hemostasis achieving patches provides a heightened level of reassurance to both the operator and patient, demonstrated through studies that cite **significant reduction in average time to deflation by 115%³ and average time to discharge by 29 – 33%.^{3,8}** Additionally, authors cite that preference from the clinical staff is given to procedures that are quicker in duration and require less manipulation of the device.¹⁷

Significant reduction in average time to deflation by 115%³ and average time to discharge by 29–33%.^{3,8}

The quicker time to hemostasis is facilitated through an increased clotting capability driven by the interaction of the Chitosan based patch independently of the normal clotting cascade. This benefit is especially apparent in patients with bleeding disorders or those who are on anticoagulant medications.^{5,6} Authors studying various clotting enhancing techniques cite significant increases in the ability to achieve initial hemostasis¹³, decreased average compression times¹, and reduced time to hemostasis.^{2,3,18-20}

In an analysis of Chitosan based patch used in conjunction with radial band technology, authors cite a 35% reduction in time to hemostasis²² and state that 69% of patients achieve hemostasis in 10 minutes or less.²¹ Fewer patient complications are also cited to be correlated with reduced compression time and strength.^{1,9-13}

Specifically in femoral artery access site procedures, authors state that **the time to hemostasis is improved by 53%¹⁹ and that 90% of patients were able to achieve hemostasis in 5 minutes or less²⁰ with the use of a Chitosan based patch** vs. standard gauze manual compression.

The incidence of infections during percutaneous coronary interventions (0.45% - 1.75%)^{23,24} is further reduced with the additional contribution of the antimicrobial properties of Chitosan.⁴ The radial artery access site provides quicker time to ambulation and authors cite little to no discomfort associated with either the band technology or the patch.²¹

In a comparison of four hemostasis achieving dressings, it was found that Celox, Quikclot, and ActCel were only 50%, 80% and 90% effective, respectfully, and that the Chitosan based patch which is used in the radial band technology was 100% effective in achieving and maintaining hemostasis.²⁵ Additionally, this patch technology can be employed across a wide range of bleeding control events including: maxillofacial surgery^{6,26-30}, hemodialysis procedures^{5,31}, penetrating trauma^{30,32,33}, emergency medical services (EMS)^{25,34,35} and combat related injuries, where Chitosan is shown to be effective against severe hemorrhage.³⁶⁻³⁹

Each year, about 1 in 25 U.S. hospital patients is diagnosed with at least one infection related to hospital care alone.⁴⁷

Estimated cost associated with hospital acquired infection is approximately \$20K.⁴⁸

**100%
effective**

Celox, Quikclot, and ActCel were only 50%, 80% and 90% effective, respectfully, and that the **Chitosan based pad which is used in the radial band technology was 100% effective** in achieving and maintaining hemostasis.²⁵

Radial vs Femoral Percutaneous Coronary Intervention (PCI) Access

The adoption rate of the transradial artery (TRA) approach as an adjunct procedure to transfemoral artery (TFA) access across the U.S. is slowly approaching the utilization rate found across Western Countries. The European Society of Cardiology recommends the use of the TRA approach as a class 1A recommendation in experienced radial operators over the TFA approach for non-ST elevated myocardial infarctions (NSTEMI).^{2,40}


In an observational study involving over 593,000 patients in the American College of Cardiology (ACC NCDR) CathPCI Registry undergoing femoral or radial procedures, authors demonstrated that **the radial approach was associated with a 67% reduction in bleeding and vascular complications** as compared to the femoral approach, without an increase in procedural failure.⁴¹

Large meta-analysis and contemporary randomized controlled trials have shown that the risk of major bleeding, all-cause mortality and major adverse cardiovascular events (MACE) are significantly reduced with use of the TRA approach as compared to TFA access.²³

Lower Cost Associated with TRA access as compared to TFA

As compared to TFA access, the TRA approach is associated with a decreased cost of hospitalization (difference of \$830) and total unadjusted procedural cost (difference of \$916).^{42,43} Additionally, same day discharge procedures are associated with cost differences as high as \$3,502.⁴³

As compared to TFA access, the TRA approach is associated with a decreased cost of hospitalization (difference of \$830) and total unadjusted procedural cost (difference of \$916).^{42,43}



References

1. Roberts JS, Niu J, Pastor-Cervantes JA. Comparison of Hemostasis Times with a Chitosan-Based Hemostatic Pad (Clo-Sur(Plus) Radial) vs Mechanical Compression (TR Band(R)) Following Transradial Access: A pilot Study. *Cardiovasc Revasc Med* 2019;20:871-4.
2. Ayyaz UI Haq M, Nazir SA, Rashid M, et al. Accelerated patent hemostasis using a procoagulant disk: a protocol designed to minimize the risk of radial artery occlusion following cardiac catheterization. *Cardiovasc Revasc Med* 2019;20:137-42.
3. Seto AH, Rollefson W, Patel MP, et al. Radial haemostasis is facilitated with a potassium ferrate haemostatic patch: the Statseal with TR Band assessment trial (STAT). *EuroIntervention* 2018;14:e1236-e42.
4. Burkatovskaya M, Tegos GP, Swietlik E, Demidova TN, A PC, Hamblin MR. Use of chitosan bandage to prevent fatal infections developing from highly contaminated wounds in mice. *Biomaterials* 2006;27:4157-64.
5. Misgav M, Kenet G, Martinowitz U. Chitosan-based Dressing for the Treatment of External/Accessible Bleedings in Children with Bleeding Tendency. *J Pediatr Hematol Oncol* 2013;00:1-3.6.
6. Kale TP, Singh AK, Kotrashetti SM, Kapoor A. Effectiveness of Hemcon dental dressing versus conventional method of haemostasis in 40 pateints on oral antiplatelet drugs. *Clinical & Basic Research* 2012;12:330-5.
7. Pappas M, Riley T, Matella T, et al. Initial Experience with the HemCon Bandage for Vascular Access Management in a Diverse Pediatric and Adult Congenital Heart Disease Patient Population. 2008: Rush Center for Congenital and Structural Heart Disease.
8. Van Meter C, Vasudevan A, Cuccerre JM, Schussler JM. Time to discharge following diagnostic coronary procedures via transradial artery approach: A comparison of Terumo band and StatSeal hemostasis. *Cardiovascular Revascularization Medicine* 2018;19:759-61.
9. Rashid M, Kwok CS, Pancholy S, et al. Radial Artery Occlusion After Transradial Interventions: A Systematic Review and Meta-Analysis. *J Am Heart Assoc* 2016;5.
10. Pancholy S, Coppola J, Patel T, Roke-Thomas M. Prevention of RAO patent hemostasis evaluation trial (PROPHET Study): A randomized coparison of traditional versus patency documented hemostasis after transradial catheterization. *Catheterization and cardiovascular interventions* 2008;72:335-40.
11. Pancholy SB, Bernat I, Bertrand OF, Patel TM. Prevention of Radial Artery Occlusion After Transradial Catheterization: The PROPHET-II Randomized Trial. *JACC Cardiovasc Interv* 2016;9:1992-9.
12. Cubero JM, Lombardo J, Pedrosa C, et al. Radial compression guided by mean artery pressure versus standard compression with a pneumatic device (RACOMAP). *Catheter Cardiovasc Interv* 2009;73:467-72.
13. Roberts JS, Niu J, Pastor-Cervantes JA. Comparison of hemostasis times with kaolin based hemostatic pad (QuikClot Radial) vs. Mechanical compression (TR Band) Following transradial access, a pilot prospective study. *J Invasive Cardiol* 2017;29:1-7.
14. Petroglou D, Didagelos M, Chalikias G, et al. Manual Versus Mechanical Compression of the Radial Artery After Transradial Coronary Angiography: The MEMORY Multicenter Randomized Trial. *JACC Cardiovasc Interv* 2018;11:1050-8.
15. Kranokpiraksa P, Pavcnik D, Kakizawa H, et al. Hemostatic efficacy of chitosan-based bandage for closure of percutaneous arterial access sites: An experimental study in heparinized sheep model. *Radiol Oncol* 2010;44:86-91.
16. Sanghvi KA, Montgomery M, Varghese V. Effect of hemostatic device on radial artery occlusion: A randomized comparison of compression devices in the radial hemostasis study. *Cardiovasc Revasc Med* 2018;19:934-8.
17. Politi L, Aprile A, Paganelli C, et al. Randomized clinical trial on short-time compression with Kaolin-filled pad: a new strategy to avoid early bleeding and subacute radial artery occlusion after percutaneous coronary intervention. *J Interv Cardiol* 2011;24:65-72.
18. Condry H, Jara C. Use of StatSeal advanced disc to decrease time to hemostasis in transradial cardiac procedures a quality improvement project. *Internation Journal of Nursing Science* 2016;6:103-7.
19. Oozawa S. A New Hemostasis Tool after Percutaneous Angioplasty: The Hemcontm Pad Hemostasis Device. *Journal of Vascular Medicine & Surgery* 2014;02.
20. Arbel J, Rozenbaum E, Reges O, et al. Usage of chitosan for femoral (USF) haemostasis after percutaneous procedures: a comparative open label study. *Euro Intervention* 2010;EPub.
21. Kang SH, Han D, Kim S, et al. Hemostasis pad combined with compression device after transradial coronary procedures: A randomized controlled trial. *PLoS One* 2017;12:e0181099.
22. Dai N, Xu DC, Hou L, Peng WH, Wei YD, Xu YW. A comparison of 2 devices for radial artery hemostasis after transradial coronary intervention. *J Cardiovasc Nurs* 2015;30:192-6.
23. Carrozza JP, Levin T. Periprocedural complications of percutaneous coronary intervention. *UpToDate*: Wolters Kluwer; 2020.
24. Hamsch JZ, Kerfeld MJ, Kirkpatrick DR, et al. Arterial Catheterization and Infection: Toll-like Receptors in Defense against Microorganisms and Therapeutic Implications. *Clinical and Translational Science* 2015;8:857-70.
25. MacIntyre AD, Quick JA, Barnes SL. Hemostatic Dressings Reduce Tourniquet Time While Maintaining Hemorrhage Control. *American Association for the Surgery of Trauma*; 2009; Pittsburgh, PA.
26. Pippi R, Santoro M, Cafolla A. The effectiveness of a new method using an extra-alveolar hemostatic agent after dental extractions in older patients on oral anticoagulation treatment: an inpatient study. *Oral Surg Oral Med Oral Pathol Oral Radiol* 2015;120:15-21.
27. Eldibany RM. Platelet rich fibrin versus Hemcon dental dressing following dental extraction in patients under anticoagulant therapy. *Tanta Dental Journal* 2014;11:75-84.
28. Azargoon H, Williams BJ, Solomon ES, Kessler HP, He J, Spears R. Assessment of hemostatic efficacy and osseous wound healing using HemCon dental dressing. *J Endod* 2011;37:807-11.
29. Malmquist JP, Clemens SC, Oien HJ, Wilson SL. Hemostasis of oral surgery wounds with the HemCon Dental Dressing. *J Oral Maxillofac Surg* 2008;66:1177-83.
30. Dailey RA, Chavez MR, Choi D. Use of a chitosan-based hemostatic dressing in dacryocystorhinostomy. *Ophthalmic Plast Reconstr Surg* 2009;25:350-3.
31. Bachtell N, Goodell T, Grunkememeier G, Jin R, Gregory K. Treatment of dialysis access puncture wound bleeding with chitosan dressings. *Dialysis & Transplantation* 2006;1-6.
32. Acheson EM, Kheirabadi BS, Deguzman R, Dick EJ, Jr., Holcomb JB. Comparison of hemorrhage control agents applied to lethal extremity arterial hemorrhages in swine. *J Trauma* 2005;59:865-74; discussion 74-5.
33. Alam HB, Chen Z, Jaskille A, et al. Application of a zeolite hemostatic agent achieves 100% survival in a lethal model of complex groin injury in Swine. *J Trauma* 2004;56:974-83.
34. Brown MA, Daya MR, Worley JA. Experience with chitosan dressings in a civilian EMS system. *J Emerg Med* 2009;37:1-7.
35. Bulger EM, Snyder D, Schoelles K, et al. An evidence-based prehospital guideline for external hemorrhage control: American College of Surgeons Committee on Trauma. *Prehosp Emerg Care* 2014;18:163-73.
36. Bennett BL, Littlejohn LF, Kheirabadi BS, et al. Management of external hemorrhage in tactical combat casualty care, chitosan-based hemostatic gauze dressings: U.S. Army Institute of Surgical Research 2014.
37. Mabry R, McManus JG. Prehospital advances in the management of severe penetrating trauma. *Crit Care Med* 2008;36:S258-66.
38. McManus JG, Wedmore I. Modern Hemostatic Agents for Hemorrhage Control - A Review and Discussion of Use in Current Combat Operations. *Business Briefing: Emergency Medicine Review* 2005.
39. Wedmore I, McManus JG, Pusateri AE, Holcomb JB. A special report on the chitosan-based hemostatic dressing: experience in current combat operations. *J Trauma* 2006;60:655-8.
40. Bernat I, Aminian A, Pancholy S, et al. Best Practices for the Prevention of Radial Artery Occlusion After Transradial Diagnostic Angiography and Intervention: An International Consensus Paper. *JACC Cardiovasc Interv* 2019;12:2235-46.
41. Rao SV, Ou F-S, Wang TY, et al. Trends in the Prevalence and Outcomes of Radial and Femoral Approaches to Percutaneous Coronary Intervention. *JACC: Cardiovascular Interventions* 2008;1:379-86.
42. Amin AP, House JA, Safley DM, et al. Costs of transradial percutaneous coronary intervention. *JACC Cardiovasc Interv* 2013;6:827-34.
43. Amin AP, Patterson M, House JA, et al. Costs Associated With Access Site and Same-Day Discharge Among Medicare Beneficiaries Undergoing Percutaneous Coronary Intervention. *JACC: Cardiovascular Interventions* 2017;10:342-51.
44. American College of Cardiology. Post PCI Sheath Removal Protocol. 2018.
45. Data on file.
46. Barnes, G.D., Lucas, E., Alexander, C., and Goldberger, S.D. National Trends in Ambulatory Oral Anticoagulant Use. *The American Journal of Medicine*. 2015. 128; 1300-1305
47. Centers for Disease Control and Prevention. CDC Winnable Battles Final Report, Healthcare-Associated Infections (HAIs). U.S. Department of Health & Human Services. Accessed on 6/24/20201 from: <https://www.cdc.gov/winnablebattles/report/HAls.html>
48. Scott RD 2nd, Culler SD, Rask KJ. Understanding the Economic Impact of Health Care-Associated Infections: A Cost Perspective Analysis. *J Infus Nurs*. 2019 Mar/Apr;42(2):61-69



Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093
 Medline United States
 1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada

1-800-396-6996 | medline.ca | canada@medline.com

Medline México

01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    