



**Vascular Access Solutions**  
Delivering excellence



# ACHIEVE OPTIMAL OUTCOMES.

You demand the best. From yourself and the companies you partner with. Products that help your hospital run more effectively and efficiently. Outcomes that meet and exceed expectations, and enhance patients' lives.

We address your needs, offering innovative, clinically equivalent products for all your diagnostic and interventional procedures, including:

- » Individually packaged, single-sterile items that help you provide meticulous care.
- » Manifold Kits that can be customized to your specifications.
- » A broad portfolio of Fluid Management Solutions designed to meet lab challenges.
- » Waste Management Kits that help create efficiencies in your perioperative environment.
- » Introducer and Mini (Micro) Kits that maximize clinical efficiency while optimizing savings.
- » Custom Procedure Packs and Trays that save time and help prevent errors.

Our team of clinical specialists will evaluate your current setup, identify innovative ways to help you improve performance, work with you to develop relevant solutions, and help get your products to you on time, every time.





A blurred background image of medical equipment, including a blue IV drip chamber in the foreground, a monitor displaying vital signs (heart rate, blood pressure, and temperature) in the background, and various medical tubes and wires. The overall color palette is dominated by blues and greens.

# MAXIMIZE VALUE. MINIMIZE WASTE.


Our product portfolio includes everything you need for fluid management in your Cath & Electrophysiology labs and Interventional Radiology. From manifold kits to enhanced custom packs, our team of specialists work with you to tailor solutions that help improve procedural efficiency at a value.

- » Custom options to satisfy your needs — We'll evaluate your current setup and identify ways to help you improve performance
- » Clinical Specialists who know your business — Technical and professional expertise to solve your challenges
- » Innovative Vascular Access Solutions — We listen to our customers and keep up with industry trends to provide Healthcare Professionals with relevant, cutting-edge product design





# CUSTOMIZE YOUR KIT.



## Table of Contents

High-quality product portfolio\* built around your clinical needs.

Standard Manifold Kits	7	Syringes	18-19
Introducer Kits	8-9	Fluid Administration Sets	20
Homeostasis & Mini Access Kits	10	Contrast Sets	20
Waste Management	11	Transducers	21
Standard Access Packs	12-13	Guidewires	22
Stopcocks	14	Safety Devices	23
Manifolds	15	Low Pressure Tubing	24
High Pressure Tubing	16	Angiographic & CT Syringes	25-26
Pressure Monitor Tubing	17	Accessories	27
Introducer Needles	17		

\*The products included in this catalog are not made with natural rubber latex.




## Standard Manifold Kits

Whether you're setting the foundation for a new cath lab or just need something quickly off the shelf, these standard kits deliver performance and value.

- » Pressure monitoring tubing with clear fixed luer lock and male/female fitting (60")
- » Waste bag with 2 -line vented spike and roller clip
- » Contrast media set with standard bore vented spike
- » Narrow barrel control syringe with rotating luer lock (8 ml)

Item No.	3-Port Manifold Right Hand	Transducer	Pkg.
VASCSHPOFF1	Off, 550 PSI	-	12 ea/cs
VASCSLPON1	On, 250 PSI, Low Torque	12", w/Wing, 1- and 3-way stopcock	12 ea/cs
VASCSLPOFF1	Off, 500 PSI, Low Torque	12", w/Wing, 1- and 3-way stopcock	12 ea/cs



## Radial Introducer Kits

- » Uniquely smooth transitions between sheath and dilator, dilator and guidewire
- » Advanced dilator snaps into place to clearly indicate when it's properly connected
- » High-performance hydrophobic-coated sheath is kink-resistant
- » Color-coded ergonomic hub
- » Reliable hemostatic seal
- » 0.018" (0.460 mm) x 45 cm guide wire with a flexible tip
- » 3-way stopcock, 21 gauge needle with protector
- » 10 ml polypropylene syringe

Item No.	Size (Fr)	Sheath Length	Dilator Length	Pkg.
DYNJINTROR4	4	5 cm	10 cm	5 ea/cs
DYNJINTROR5	5	5 cm	10 cm	5 ea/cs
DYNJINTROR6	6	5 cm	10 cm	5 ea/cs

### Did you know?

Medline is ISO 13485 certified—a distinction that's granted only after rigorous analysis of a company's quality systems.



## Femoral/Subclavian Introducer Kits

Used for the introduction in the femoral or subclavian artery for vein or vascular procedures.

- » Uniquely smooth transitions between sheath and dilator, dilator and guidewire
- » Advanced dilator snaps into place when it's properly connected
- » High-performance hydrophobic-coated sheath is kink-resistant
- » Color-coded ergonomic hub
- » Reliable hemostatic seal
- » 0.038" (0.97 mm) x 50 cm guidewire with flexible straight and "J" tips
- » 3-way stopcock, 18-gauge needle with protector
- » 10 ml polypropylene syringe

### Standard Kit

Item No.	Size (Fr)	Sheath Length	Dilator Length	Pkg.
DYNJINTROS4	4	10 cm	15 cm	5 ea/cs
DYNJINTROS5	5	10 cm	15 cm	5 ea/cs
DYNJINTROS6	6	10 cm	15 cm	5 ea/cs
DYNJINTROS7	7	10 cm	15 cm	5 ea/cs
DYNJINTROS8	8	10 cm	15 cm	5 ea/cs
DYNJINTROS9	9	10 cm	15 cm	5 ea/cs
DYNJINTROS10	10	10 cm	15 cm	5 ea/cs
DYNJINTROS11	11	10 cm	15 cm	5 ea/cs

### Valved Kit

Item No.	Size (Fr)	Sheath Length	Dilator Length	Pkg.
DYNJINTROL4	4	23 cm	28 cm	5 ea/cs
DYNJINTROL5	5	23 cm	28 cm	5 ea/cs
DYNJINTROL6	6	23 cm	28 cm	5 ea/cs
DYNJINTROL7	7	23 cm	28 cm	5 ea/cs
DYNJINTROL8	8	23 cm	28 cm	5 ea/cs
DYNJINTROL9	9	23 cm	28 cm	5 ea/cs
DYNJINTROL10	10	23 cm	28 cm	5 ea/cs
DYNJINTROL11	11	23 cm	28 cm	5 ea/cs
DYNJINTROL12	12	23 cm	28 cm	5 ea/cs

### Peel Away Kit

Item No.	Size (Fr)	Sheath Length	Dilator Length	Pkg.
DYNJINTROP6	6	16 cm	19 cm	5 ea/cs
DYNJINTROP7	7	16 cm	19 cm	5 ea/cs
DYNJINTROP75	7.5	16 cm	19 cm	5 ea/cs
DYNJINTROP8	8	16 cm	19 cm	5 ea/cs
DYNJINTROP85	8.5	16 cm	19 cm	5 ea/cs
DYNJINTROP9	9	16 cm	19 cm	5 ea/cs
DYNJINTROP95	9.5	16 cm	19 cm	5 ea/cs
DYNJINTROP10	10	16 cm	19 cm	5 ea/cs
DYNJINTROP105	10.5	16 cm	19 cm	5 ea/cs
DYNJINTROP11	11	16 cm	19 cm	5 ea/cs
DYNJINTROP115	11.5	16 cm	19 cm	5 ea/cs
DYNJINTROP12	12	16 cm	19 cm	5 ea/cs
DYNJINTROP125	12.5	16 cm	19 cm	5 ea/cs





## Mini Access Kits

Enable successful placement during vein or vascular procedures.

- » 0.18" x 45 cm guidewire dia/length
- » 10 cm introducer sheath length

### Stainless Steel Guidewire, Regular Style

Item No.	Sheath Size (Fr)	Introducer Needle	Pkg.
DYNJMAST01	4	21 g x 2.75", No base	5 ea/cs
DYNJMAST02	5	21 g x 2.75", No base	5 ea/cs

### Stainless Steel Guidewire, Stiff Style

Item No.	Sheath Size (Fr)	Introducer Needle	Pkg.
DYNJMAST07	4	21 g x 2.75", Echogenic	5 ea/cs
DYNJMAST08	5	21 g x 2.75", Echogenic	5 ea/cs

### Nitinol/Tungsten Guidewire, Regular Style

Item No.	Sheath Size (Fr)	Introducer Needle	Pkg.
DYNJMAST03A	4	21 g x 2.75", No base	5 ea/cs
DYNJMAST04	5	21 g x 2.75", No base	5 ea/cs

### Nitinol/Tungsten Guidewire, Stiff Style

Item No.	Sheath Size (Fr)	Introducer Needle	Pkg.
DYNJMAST09	4	21 g x 2.75", Echogenic	5 ea/cs
DYNJMAST10	5	21 g x 2.75", Echogenic	5 ea/cs



## Hemostatic Device Kit

- » Allows for a secure hemostasis seal with the capability to use interventional componentry
- » Simple, single-handed maneuverability Y-connector
- » Insertion tool 20 g x 10.9 cm
- » Torque device ABS material accepts guidewires with 0.014 to 0.022" diameters

Item No.	Pkg.	Nonsterile No.
DYNJHEMOKIT	25 ea/cs	91763

### Did you know?

We've been serving facilities like yours since 1966.



## Closed Waste Management

Designed to collect patient or operating waste to minimize staff exposure to biohazard waste and bloodborne pathogens.

- » Minimize splashing and limit potential for exposure to medical waste
- » Customizable to your needs
- » Available in single sterile or as part of a custom kit
- » Can be used off a manifold or as a closed system
- » Includes a drape clip for security during procedures

### 1-Line Clear Waste Bag

Item No.	Fitting	Reflux Valve	3-Way Stopcock	Pkg.	Nonsterile No.
DYNJWB DIR	Male	–	–	20 ea/cs	89345

### 2-Line Clear Waste Bag, Clear Vented Spike

Item No.	Fitting	Reflux Valve	3-Way Stopcock	Pkg.	Nonsterile No.
DYNJWBSSRVA	Male	●	–	20 ea/cs	83494/ 94412
DYNJWBSSRV	Female	●	–	20 ea/cs	83494
DYNJWBSSA	Male	–	–	20 ea/cs	48900/94412
DYNJWBSS	Female	–	–	20 ea/cs	48900

### 2-Line Clear Waste Bag, Drip Chamber, Clear Filtered


Item No.	Fitting	Reflux Valve	3-Way Stopcock	Pkg.	Nonsterile No.
DYNJWB DCA	Male	–	–	20 ea/cs	98975/ 94412
DYNJWBDC	Female	–	–	20 ea/cs	98975

### 3-Line Clear Waste Bag, Clear Vented Spike, Green Vented Spike

Item No.	Fitting	Reflux Valve	3-Way Stopcock	Pkg.	Nonsterile No.
DYNJWBDSRVA	Male	●	●	20 ea/cs	100435/ 94412
DYNJWBDSRV	Female	●	●	20 ea/cs	100435
DYNJWBDSA	Male	–	●	20 ea/cs	83493/ 94412
DYNJWBDS	Female	–	●	20 ea/cs	83493

### 3-Line Clear Waste Bag, Clear Vented Spike, Clear Port, 2 Check Valves

Item No.	Fitting	Reflux Valve	3-Way Stopcock	Pkg.	Nonsterile No.
DYNJWB CRVA	Male	●	●	20 ea/cs	92110/ 94412
DYNJWB CRV	Female	●	●	20 ea/cs	92110



**DRAPE CLIP**



- » Standard sterile packs are stocked and ready to order at any time
- » Any item needed for your cardiology or interventional radiology procedures can be added to a custom pack by working with your Medline representative
- » The high volume of standard packs helps keep costs low

Item No.	Description	Pkg.
DYNJCATH1	Catheter Lab Pack 1	25/cs

Qty (ea)	Description	Qty (ea)	Description	Qty (ea)	Description
1	Band Bag, 25" Deep, 50" Diameter	50	Gauze, 4" x 4", 12-ply	1	PVP Solution, 2-oz., Paint
1	Band Bag, 36" x 28" with Tape	1	Vinyl Exam Glove, L	1	PVP Solution, 2-oz., Scrub
2	Bowl, 32-oz., Pigment-Free	2	Poly Reinforced Gown, XL	4	Winged Sponge, 1.75" x 3.5" x 1"
1	Bowl, Guide Wire, with Tabs	1	Vascular Access Label Sheet	5	10 ml Syringe, Luer Lock
4	Towel Clamp	1	Marker, Skin/Utility Dual	1	20 ml Syringe, Luer Lock
1	Set-Up Cover, 54" x 80"	1	Needle, 18-Gauge x 1.5"	2	3 ml Syringe, Luer Lock
1	Medicine Cup, 2-oz., Clear	1	Needle, 22-Gauge x 1.5"	1	Absorbent Towel, 17" x 20", White
1	Table Cover, 50" x 90"	1	Needle, 25-Gauge x .63"	6	OR Towel, Blue
1	Fem Angio Drape with Pouch, 80" x 135"	1	Arterial Needle, 18-Gauge x 2.75", No Wings	1	2-Compartment Tray, Pigment Free
1	Curved Mosquito Forceps	1	Safety Scalpel, #11, Stainless Steel	1	CSR Wrap, 25" x 25", White





## Standard Vascular Access Pack 2

Item No.	Description	Pkg.
DYNJSCL300	Catheter Lab Pack 2	4/cs

### Bill of Materials

Qty (ea)	Description	Qty (ea)	Description	Qty (ea)	Description
1	Male-to-Male Adaptor	1	Curved Mosquito Forceps	1	Safety Scalpel, #11 Stainless Steel
1	Waste Bag	40	Gauze, 4" x 4", 12-ply	1	OR Scissors
1	Band Bag, 25" Deep, 50" Diameter	1	Gown, Poly Reinforced, with Absorbent Towel, XL	1	IV Tubing
1	Band Bag, 36" x 28", with Tape	2	Gown, Poly Reinforced, with Breathable Sleeve, XL	1	10 ml Syringe, Luer Lock, Blue Plunger
1	Bowl, 32-oz., Pigment Free	1	Guide Wire	1	10 ml Syringe, Luer Lock, Green Plunger
1	Bowl, Guide Wire, with Tabs	1	Webster Needle Holder	1	10 ml Syringe, Luer Lock, Red Plunger
1	Chloraprep® Applicator, 10.5 ml	1	Labels	1	10 ml Syringe, Luer Lock, White Plunger
4	Towel Clamp	1	Manifold, 3-Port, 500 PSI	1	10 ml Syringe, Luer Lock, Yellow Plunger
1	Setup Cover, 54" x 80"	1	Marker, Skin/Utility Dual	1	12 ml Control Syringe, Luer Lock
1	Medicine Cup, 2-oz., Blue	1	Needle, 18-Gauge x 1.5"	1	Pressure Control Syringe, 10 ml, Luer Lock
1	Medicine Cup, 2-oz., Clear	1	Needle, 22-Gauge x 1.5"	1	Pressure Tubing
1	Medicine Cup, 2-oz., Green	1	Needle, 25-Gauge x 1.5"	1	Pressure Tubing Monitor
1	Medicine Cup, 2-oz., Red	1	Needle Nest	6	Blotting Towel
1	Medicine Cup, 2-oz., Yellow	1	Hypodermic Needle, 18-Gauge x 2.75", No Base	1	Table Cover, 50" x 90"
1	Fem Angio Drape, with Pouch, 77" x 135"				

## LOW PRESSURE



DYNJSC302



DYNJSC301



DYNJSC101

## MEDIUM PRESSURE



DYNJSC303



DYNJSC304

## HIGH PRESSURE



DYNJSC305HP



DYNJSC306HP

## Stopcocks

An externally operated valve regulating the flow of a liquid or gas.

- » Can be attached to any configuration inside a pack for custom procedures
- » Efficient on/off valve can be used in a variety of settings

### One-Way Body Style, 100 ml/minute

Item No.	Handle Position	Pressure	PSI	Male Fitting	Pkg.	Nonsterile No.
DYNJSC101	On	Low	200	Rotating Luer Lock	50 ea/cs	81205

### Three-Way Body Style, 100 ml/minute

Item No.	Handle Position	Pressure	PSI	Male Fitting	Pkg.	Nonsterile No.
DYNJSC301	Off	Low	45	Rotating Collar	20 ea/cs	90126
DYNJSC302	On	Low	45	Rotating Collar	20 ea/cs	91627
DYNJSC307	Off	Low	45	Fixed Luer Lock	20 ea/cs	07065
DYNJSC303	Off	Med	500	Rotating Luer Lock	20 ea/cs	91628
DYNJSC304	On	Med	500	Rotating Luer Lock	20 ea/cs	90127
DYNJSC305HP	Off	High	1,050	Rotating Luer Lock	20 ea/cs	91629
DYNJSC306HP	On	High	1,050	Rotating Adaptor	20 ea/cs	90128


### Four-Way Body Style, Hi-Flow, 200 ml/minute

Item No.	Handle Position	Pressure	PSI	Male Fitting	Pkg.	Nonsterile No.
DYNJSC401	Off	Low	45	Rotating Luer Lock	50 ea/cs	81203

### Did you know?

We offer over 35,000 components to meet every healthcare need.





## Manifolds

Manifolds are the primary vehicle used to provide access to saline, contrast, pressure monitoring and waste management systems.

- » Easy to use, smooth turning handles
- » Clearly identified On and Off options
- » Requires minimal stopcock manipulations
- » Can be placed inside a custom kit

### ON/OFF HANDLE

On



Off



### Standard Handle

Item No.	No. of Ports	Orientation	On/Off Handle	PSI	Fitting	Pkg.	Nonsterile No.
DYNJMAN01	2	Right Hand	On	500	Rotating Male Luer Lock-Female	10 ea/cs	91622
DYNJMAN02	2	Right Hand	On	500	Rotating Male Luer Lock-Female	10 ea/cs	91623
DYNJMAN03	2	Right Hand	Off	500	Rotating Male Luer Lock-Female	10 ea/cs	91624
DYNJMAN04	3	Right hand	Off	500	Rotating Male Luer Lock-Female	10 ea/cs	91625
DYNJMAN05	4	Right Hand	On	500	Rotating Male Luer Lock-Female	10 ea/cs	91626
DYNJMAN06	4	Right Hand	Off	500	Rotating Male Luer Lock-Female	10 ea/cs	90125
DYNJMAN32	4	360°	Off	29	Female/Female	25 ea/cs	112079

## High Pressure Contrast Injection Tubing

Contrast injection lines are used to connect components together, pass fluids or medication through the manifold, or as extensions between other medical elements.

- » Sturdy PVC material
- » Clear tubing has secured connections to mitigate against air in the line and is rated to 1,050 PSI
- » Braided polyurethane design allows for flexibility of tube while maintaining consistent performance, is flexible to accommodate the preferences of your physicians, and withstands 1,200 PSI.

### High Pressure, Clear Tubing – ID .071" (1.8 mm), OD .019" (4.75 mm)

Item No.	Length	Fitting	Pkg.	Nonsterile No.
DYNJTBGHP10F	10" (25.40 cm)	Male Fixed Luer Lock – Female	50 ea/cs	93058
DYNJTBGHP10R	10" (25.40 cm)	Male Rotating Luer Lock – Female	50 ea/cs	93059
DYNJTBGHP20F	20" (50.80 cm)	Male Fixed Luer Lock – Female	50 ea/cs	93060
DYNJTBGHP20R	20" (50.80 cm)	Male Rotating Luer Lock – Female	50 ea/cs	93061
DYNJTBGHP30F	30" (76.20 cm)	Male Fixed Luer Lock – Female	50 ea/cs	93062
DYNJTBGHP30R	30" (76.20 cm)	Male Rotating Luer Lock – Female	50 ea/cs	93063
DYNJTBGHP48F	48" (121.92 cm)	Male Fixed Luer Lock – Female	50 ea/cs	93064
DYNJTBGHP48R	48" (121.92 cm)	Male Rotating Luer Lock – Female	50 ea/cs	93065
DYNJTBGHP60F	60" (152.40 cm)	Male Fixed Luer Lock – Female	50 ea/cs	93066
DYNJTBGHP60R	60" (152.40 cm)	Male Rotating Luer Lock – Female	50 ea/cs	93067
DYNJTBGHP72F	72" (182.88 cm)	Male Fixed Luer Lock – Female	50 ea/cs	93068
DYNJTBGHP72R	72" (182.88 cm)	Male Rotating Luer Lock – Female	50 ea/cs	93069
DYNJTBGHP84F	84" (213.36 cm)	Male Fixed Luer Lock – Female	50 ea/cs	93070
DYNJTBGHP84R	84" (213.36 cm)	Male Rotating Luer Lock – Female	50 ea/cs	93071
DYNJTBGHP96F	96" (243.84 cm)	Male Fixed Luer Lock – Female	50 ea/cs	93072
DYNJTBGHP96R	96" (243.84 cm)	Male Rotating Luer Lock – Female	50 ea/cs	93073

### High Pressure, Braided Tubing – ID .071" (1.8 mm), OD .014" (3.65 mm)

Item No.	Length	Fitting	Pkg.	Nonsterile No.
DYNJBRT10F	10" (25.40 cm)	Male Fixed Luer Lock – Female	50 ea/cs	90511
DYNJBRT10R	10" (25.40 cm)	Male Rotating Luer Lock – Female	50 ea/cs	90512
DYNJBRT20F	20" (50.80 cm)	Male Fixed Luer Lock – Female	50 ea/cs	90513
DYNJBRT20R	20" (50.80 cm)	Male Rotating Luer Lock – Female	50 ea/cs	90514
DYNJBRT30F	30" (76.20 cm)	Male Fixed Luer Lock – Female	50 ea/cs	90515
DYNJBRT30R	30" (76.20 cm)	Male Rotating Luer Lock – Female	50 ea/cs	90516
DYNJBRT48F	48" (121.92 cm)	Male Fixed Luer Lock – Female	50 ea/cs	90517
DYNJBRT48R	48" (121.92 cm)	Male Rotating Luer Lock – Female	50 ea/cs	90518
DYNJBRT60F	60" (152.40 cm)	Male Fixed Luer Lock – Female	50 ea/cs	90519
DYNJBRT60R	60" (152.40 cm)	Male Rotating Luer Lock – Female	50 ea/cs	90520
DYNJBRT72F	72" (182.88 cm)	Male Fixed Luer Lock – Female	50 ea/cs	90521
DYNJBRT72R	72" (182.88 cm)	Male Rotating Luer Lock – Female	50 ea/cs	90522

#### LUER LOCK STYLES

##### Rotating Luer Lock



##### Fixed Luer Lock



Detail of braided tubing.

## Pressure Monitoring Tubing

Used to transmit internal blood pressure from the patient to monitoring systems.

- » Non-expanding PVC material
- » Ensures high-quality wave forms with improved clarity

Item No.	Length	Fitting	Pkg.	Nonsterile No.
DYNJPMTBG6MF	6" (15.24 cm)	Male Fixed Luer Lock – Female	25 ea/cs	120910
DYNJPMTBG12MM	12" (30.48 cm)	Male – Male Fixed Luer Locks	25 ea/cs	129011
DYNJPMTBG12MF	12" (30.48 cm)	Male Fixed Luer Lock – Female	25 ea/cs	129012
DYNJPMTBG24MM	24" (60.96 cm)	Male – Male Fixed Luer Locks	25 ea/cs	129013
DYNJPMTBG24MF	24" (60.96 cm)	Male Fixed Luer Lock – Female	25 ea/cs	129014
DYNJPMTBG36MM	36" (91.44 cm)	Male – Male Fixed Luer Locks	25 ea/cs	120915
DYNJPMTBG36MF	36" (91.44 cm)	Male Fixed Luer Lock – Female	25 ea/cs	120916
DYNJPMTBG48MM	48" (121.92 cm)	Male – Male Fixed Luer Locks	25 ea/cs	120917
DYNJPMTBG48MF	48" (121.92 cm)	Male Fixed Luer Lock – Female	25 ea/cs	120918
DYNJPMTBG60MM	60" (152.40 cm)	Male – Male Fixed Luer Locks	25 ea/cs	120919
DYNJPMTBG60MF	60" (152.40 cm)	Male Fixed Luer Lock – Female	25 ea/cs	120920
DYNJPMTBG72MM	72" (182.88 cm)	Male – Male Fixed Luer Locks	25 ea/cs	120921
DYNJPMTBG72MF	72" (182.88 cm)	Male Fixed Luer Lock – Female	25 ea/cs	120922
DYNJPMTBG84MF	84" (213.36 cm)	Male Fixed Luer Lock – Female	25 ea/cs	120924
DYNJPMTBG96MF	96" (243.84 cm)	Male Fixed Luer Lock – Female	25 ea/cs	120926

*Hemodynamic cables, transducers and accessories available on page 21*

### Did you know?

We control the entire supply chain, so you win on quality, service and price.

## Angiography Introducer Needles

Used to gain access to the veins or arteries.

- » Needle points have a non-coring "B" arterial bevel
- » Single wall arterial percutaneous puncture
- » Available with a window hub and winged base plate
- » 2.75" long

Item No.	Needle Gauge	Wings/Base	Echogenic	Color	Pkg.	Nonsterile No.
DYNJNDL18	18	–	–	Pink	25 ea/cs	83499
DYNJNDL18WING	18	●	–	Pink	25 ea/cs	84880
DYNJNDL19	19	–	–	Brown	25 ea/cs	84881
DYNJNDL20	20	–	–	Yellow	25 ea/cs	89331
DYNJNDL21	21	–	–	Green	25 ea/cs	89332
DYNJNDL21E	21	–	●	Lt. Green	25 ea/cs	103763
DYNJNDL22	22	–	–	White	25 ea/cs	89346





## Angiography Control Syringes

A more secure and controlled vehicle used for aspiration, contrast and/or medication.

- » A variety of configurations for your team's clinical preference
- » Easy visualization of air bubbles with clear polycarbonate design
- » Requires less force and provides smoother injections and aspirations with improved barrel design

### SYRINGE TIP STYLES

Fixed



Rotating



Item No.	Volume (ml)	Male Fitting	Plunger Style	Narrow Barrel	Pkg.	Nonsterile No.
DYNJSYR8CWR	8	Rotating	Thumb Ring	●	10 ea/cs	89369
DYNJSYR8C	8	Fixed	Thumb Ring	●	10 ea/cs	89371
DYNJSYR10CWR	10	Rotating	Thumb Ring	—	10 ea/cs	81015
DYNJSYR10CRP	10	Rotating	Palm Pad	—	10 ea/cs	96723
DYNJSYR10CP	10	Fixed	Palm Pad	—	10 ea/cs	96725
DYNJSYR10C	10	Fixed	Thumb Ring	—	10 ea/cs	88618
DYNJSYR10CB	10	Fixed	Thumb Ring	●	10 ea/cs	120185
DYNJSYR10CWRB	10	Rotating	Thumb Ring	●	10 ea/cs	120184
DYNJSYR12C	12	Fixed	Thumb Ring	—	10 ea/cs	88619
DYNJSYR12CWR	12	Rotating	Thumb Ring	—	10 ea/cs	81016
DYNJSYR12CRP	12	Rotating	Palm Pad	—	10 ea/cs	96724
DYNJSYR12CP	12	Fixed	Palm Pad	—	10 ea/cs	96726

## Colored Syringes

Used for aspiration and injection of contrast, saline or medication.

### Customizable

- » Customize with your specific label preference
- » Choose from a variety of sizes and barrel colors to meet your specific needs
- » Available inside angio kits as well as vascular access packs
- » Additional colors and sizes available for use inside custom packs

### Durable

- » Polycarbonate withstands the high pressure ratings in the vascular access lab
- » Individually packaged, sterile
- » Additional colors available for use inside custom kits

### 1 ml Volume

Item No.	Color	Pkg.	Nonsterile No.
DYNJSYRPC1R	Red	10 ea/cs	83456
DYNJSYRPC01W	White	10 ea/cs	83453

### 3 ml Volume

Item No.	Color	Pkg.	Nonsterile No.
DYNJSYRPC3B	Blue	10 ea/cs	83551
DYNJSYRPC3G	Green	50 ea/cs	92189
DYNJSYRPC3R	Red	10 ea/cs	83553
DYNJSYRPC03W	White	10 ea/cs	83550
DYNJSYRPC3Y	Yellow	10 ea/cs	83552

### 6 ml Volume

Item No.	Status	Pkg.	Nonsterile No.
DYNJSYRPC6B	Blue	10 ea/cs	83557
DYNJSYRPC6R	Red	10 ea/cs	83558
DYNJSYRPC6W	White	10 ea/cs	83556
DYNJSYRPC6Y	Yellow	50 ea/cs	92191

### 10 ml Volume

Item No.	Color	Pkg.	Nonsterile No.
DYNJSYRPC10B	Blue	10 ea/cs	103548
DYNJSYRPC10G	Green	10 ea/cs	103551
DYNJSYRPC10P	Purple	10 ea/cs	103552
DYNJSYRPC10R	Red	10 ea/cs	103550
DYNJSYRPC10W	White	10 ea/cs	103546
DYNJSYRPC10Y	Yellow	10 ea/cs	103549



### 20 ml Volume

Item No.	Color	Pkg.	Nonsterile No.
DYNJSYRPC20B	Blue	10 ea/cs	83592
DYNJSYRPC20R	Red	10 ea/cs	83594
DYNJSYRPC20W	White	10 ea/cs	83591
DYNJSYRPC20BY	Yellow	10 ea/cs	83593

### 30 ml Volume

Item No.	Color	Pkg.	Nonsterile No.
DYNJSYRPC30B	Blue	10 ea/cs	83598
DYNJSYRPC30P	Purple	10 ea/cs	83559
DYNJSYRPC30W	White	10 ea/cs	83597

## Standard Vascular Access Label Kit

We have custom capabilities to design the labels you need.

- » In color or black and white
- » Quantities that suit your practice
- » Verbiage that aligns with your procedure
- » Ultra-fine point marker included

Item No.	Pkg.	Nonsterile No.
DYNJCATHLBL	12 ea/cs	95806 Label, 44035 Marker

Label	Label Count
Adenosine	1
Blank	10
Contrast	4
Heparin	2
Heparin/Saline	4
Lidocain	2
Nitroglycerin	2
Saline	4
Verapamil	1



## Fluid Administration Sets

Administer saline, medication or dye in operating rooms and labs.

- » For the majority of your clinical needs
- » Drip chamber assists in mitigating against air passed downstream
- » Venting assists in regulating the flow of fluid during a procedure
- » Includes a regulator that acts as an on/off valve

### Clear, Vented

Item No.	Length	Filter	Bore	Pkg.	Nonsterile No.
DYNJVNFD72	72" (182.88 cm)	–	Std.	25 ea/cs	92528
DYNJVFDC72	72" (182.88 cm)	●	Std.	25 ea/cs	92527

### Clear, Not Vented

Item No.	Length	Filter	Bore	Pkg.	Nonsterile No.
DYNJNVFD72	72" (182.88 cm)	–	Std.	25 ea/cs	92526
DYNJLBD72	72" (182.88 cm)	–	Large	25 ea/cs	35004
DYNJLBVS108	108" (274.32 cm)	●	Std.	25 ea/cs	92530

## Contrast Sets

Designed specifically to administer contrast or dye.

- » 72" (182.88 cm) long
- » Configure your contrast/dye container to fluid management device
- » Fittings are compatible with manifolds, stopcocks or syringes with a male luer lock

Item No.	Gauge	Check Valve w/Extension	Pkg.	Nonsterile No.
DYNJSPKVST72	Standard	–	25 ea/cs	92521
DYNJLBSPK72	Large	–	25 ea/cs	35003
DYNJTUMSYR	Large	●	10 ea/cs	–

## Short Spike

A required accessory for contrast management systems.

- » Single sterile; new spike assembly recommended for each bottle of contrast
- » Vented

Item No.	Flow Control	Pkg.	Nonsterile No.
DYNJSSPIKE	Pinch Clamp/Vented Spike	25 ea/cs	90509



## Disposable Pressure Transducers

Used to measure accurate and real time internal pressures.

- » Compact design
- » Clear, visible fluid path makes it easy to flush and de-bubble the line
- » Single-hand operational flush options
- » Cables available to ensure proper connections

Item No.	Cable Length	3-Way Stopcock	1-Way Stopcock	Flush Device	Fitting	Pkg.	Nonsterile No.
DYNJTRANSWG	12"	–	–	–	Male-Male	20 ea/cs	95103
DYNJTRANSMF	12"	–	–	–	Male-Female	20 ea/cs	95104
DYNJTRANSF	12"	●	–	●	Male-Female	20 ea/cs	95105
DYNJTRANS3	12"	●	●	–	Male-Female	20 ea/cs	99049
DYNJTRANS1	12"	–	●	–	Male-Male	20 ea/cs	114423
DYNJTRANS48	48"	–	●	–	Male-Female	20 ea/cs	103013
DYNJTRANS48	48"	●	–	–	Male-Male	20 ea/cs	98136

## Transducer Accessories

Economical program for cables, mounts and c-clamps.

- » We can also outfit your facility with the cables and accessories your monitors or hemodynamic machines require
- » Clamps and plates for secure attachment are provided with the cables upon conversion.
- » Reliable precision
- » Customizable configurations

Item No.	Description	Pkg.	Nonsterile No.
DYNJTRPLATE1	Plate	1 bx/cs, 1 ea/cs	97805
DYNJCLAMP1	Clamp	5 ea/cs	102128
DYNJCBLPHILH	Philips Cable	1 bx/ea	–
DYNJCBLSIEMH	Siemens Cable	1 bx/ea	–
DYNJCBLTRAMH	GE Cable	1 bx/ea	–



DYNJCLAMP1



DYNJTRPLATE1



## Diagnostic Guide Wires

Guide wires are used to facilitate catheter placement in the vascular system.

- » A complete line of diagnostic guide wires to ensure the perfect fit for all your procedural needs
- » Industry-leading pre-coated design ensures uniformity, lubricity, adhesion for reduced friction and atraumatic results
- » Built for smooth insertions and atraumatic design that performs

### Fixed Core

Item No.	Tip	Diameter	Length	Coating	Straightenable	Pkg.	Nonsterile No.
DYNJGWIRE01	3 mm J	0.035"	150 cm	PTFE	●	10 ea/cs, 1 bx/cs	13060
DYNJGWIRE07	3 mm J	0.035"	175 cm	PTFE	●	10 ea/cs, 1 bx/cs	22610
DYNJGWIRE03	3 mm J	0.035"	180 cm	PTFE	●	10 ea/cs, 1 bx/cs	88920
DYNJGWIRE20	3 mm J	0.035"	260 cm	PTFE	●	10 ea/cs, 1 bx/cs	102980
DYNJGWIRE06	3 mm J	0.038"	150 cm	PTFE	●	10 ea/cs, 1 bx/cs	17839

### Fixed Core, Firm Body

Item No.	Tip	Diameter	Length	Coating	Straightenable	Pkg.	Nonsterile No.
DYNJGWIRE05	3 mm J	0.035"	150 cm	PTFE	●	10 ea/cs, 1 bx/cs	22604
DYNJGWIRE08	3 mm J	0.038"	150 cm	PTFE	●	10 ea/cs, 1 bx/cs	22605

### Fixed Core, Double-ended

Item No.	Tip	Diameter	Length	Coating	Straightenable	Pkg.	Nonsterile No.
DYNJGWIRE16	Double Distal	0.035"	70 cm	None	–	10 ea/cs, 1 bx/cs	31038
DYNJGWIRE10	Double Distal	0.035"	150 cm	PTFE	●	10 ea/cs, 1 bx/cs	29307
DYNJGWIRE11	3 mm J	0.035"	150 cm	PTFE	●	10 ea/cs, 1 bx/cs	49744
DYNJGWIRE13	Double Distal	0.035"	180 cm	PTFE	●	10 ea/cs, 1 bx/cs	88923

### Fixed Core, Heparin Bond, Double-ended

Item No.	Tip	Diameter	Length	Coating	Straightenable	Pkg.	Nonsterile No.
DYNJGWIRE14	Bentson	0.038"	150 cm	PTFE	–	10 ea/cs, 1 bx/cs	88924

### Fixed Core, Heparin Bond

Item No.	Tip	Diameter	Length	Coating	Straightenable	Pkg.	Nonsterile No.
DYNJGWIRE04	3 mm J	0.035"	150 cm	PTFE	●	10 ea/cs, 1 bx/cs	22603
DYNJGWIRE09	Bentson	0.035"	150 cm	PTFE	–	10 ea/cs, 1 bx/cs	88925

### Fixed Core, Heavy Duty

Item No.	Tip	Diameter	Length	Coating	Straightenable	Pkg.	Nonsterile No.
DYNJGWIRE02	3 mm J	0.035"	150 cm	PTFE	●	10 ea/cs, 1 bx/cs	32736
DYNJGWIRE12	3 mm J	0.038"	180 cm	PTFE	●	10 ea/cs, 1 bx/cs	88921


### Movable Core

Item No.	Tip	Diameter	Length	Coating	Straightenable	Pkg.	Nonsterile No.
DYNJGWIRE17	3 mm J	0.038"	150 cm	PTFE	–	10 ea/cs, 1 bx/cs	13129
DYNJGWIRE19	3 mm J	0.038"	150 cm	PTFE	●	10 ea/cs, 1 bx/cs	102980

### Movable Core, Heparin Bond

Item No.	Tip	Diameter	Length	Coating	Straightenable	Pkg.	Nonsterile No.
DYNJGWIRE18	3 mm J	0.035"	150 cm	PTFE	–	10 ea/cs, 1 bx/cs	22609

### GUIDE WIRE DETAIL



Pre-coated design ensures uniform coating and reduced friction.



## Guide Wire Bowls

Provide a secure location while preparing for a procedure.

- » Easily graspable while wearing surgical gloves
- » Five tabs ensure that the wires stay in place

Item No.	Diameter	Fluid Capacity	Pkg.	Nonsterile No.
DYNJGUIDEWIRE	8.5"	2,500 ml	28 ea/cs	96360
DYNJGUIDEB2	11"	5,000 ml	10 ea/cs	119542

## Safety Devices

- » Safety splash self-contained fluid management system eliminates the need to transport medical waste in liquid form
- » Sterile foam needle holder offers temporary placement of needles

### Safety Splash 500 ml Fluid Waste Bowl with Pad, Gold

Item No.	Pkg.	Nonsterile No.
DYNJSPLASH1	24 ea/cs	114270

### Foam/Adhesive Needle Holder

Item No.	Pkg.	Nonsterile No.
DYNJNDLNEST	50 ea/cs	35862

## Fluid Transfer Devices

Medline's fluid decanters help ensure an aseptic transfer or removal of fluids from containers.

- » Single use
- » Functional in all operating room environments

Item No.	Description	Pkg.	Nonsterile No.
DYNJDEC03	Decanter 3" Transfer	50 ea/cs	47578
DYNJDEC06	Decanter 6" Vial	50 ea/cs	83480
DYNJDEC09	Decanter 9" Bag	50 ea/cs	92101







## Low Pressure Tubing

- » 1,200 PSI rated
- » Offers greater flexibility
- » Less expensive than braided tubing
- » Clear design offers better visual bubble detection
- » Coiled packaging maximizes storage space
- » Y-Lines are 330 PSI rated

### Did you know?

We have 7 pack assembly sites worldwide and thousands of employees making over 125K surgical packs each day.

### Couer Coiled Extension

Item No.	Fits	Includes	Length	Pkg.
VINC4050060	Mallinkrodt Pre-Filled Syringes	60" Coiled Extension Line	60"	100 ea/cs

### Couer Y-Lines

Item No.	Fits	Includes	Length	Pkg.
VINC4052155	Optivantage	60" Straight Y-Line with 2 Check Valves	60"	100 ea/cs
VINC4052604	Optivantage	60" Straight Y-Line with 1 Check Valve	60"	100 ea/cs
VINC4052962	Mallinkrodt Injectors	96" Straight Y-Line with 2 Check Valves	96"	100 ea/cs



## Couer CT Syringes

- » Direct Connect tip design allows for a quick, leak-free connection of transfer sets, spikes and extension lines, while reducing contrast waste
- » Extended length plunger provides trouble-free syringe filling

Item No.	Fits	Includes	Size	Pkg.
VINC8530125	Mallinckrodt Optivantage and Liebel-Flarshiem	Fill Tube	125 ml	50 ea/cs
VINC8530126	Basin With Tabs, 12"	60" Low Pressure Extension Line and Spike	125 ml	40 ea/cs
VINC8530127	Mallinckrodt Optivantage and Liebel-Flarshiem	Fill Tube and 60"	125 ml	50 ea/cs
VINC8530211	Medrad MCT, MCT Plus, Liebel-Flarshiem CT9000 Series Injectors	Low Pressure Extension Line	200 ml	50 ea/cs



## Couer Angiographic Syringes

Item No.	Fits	Includes	Size	Pkg.
VINC8530155**	Medrad Mark V, Mark V Plus, and Provis Injectors	Fill Tube	150 ml	50 ea/cs
VINC8530210	Medrad Mark V, Mark V Plus, and Provis Injectors	Fill Tube	200 ml	50 ea/cs
VINC8536	Medrad Mark V, Mark V Plus, and Provis Injectors	Fill Tube and 60"	150 ml	50 ea/cs

\*\*Adaptor/Pressure Jacket required for use on Illumena.

## Couer MRI Syringe

Item No.	Fits	Includes	Size	Pkg.
VINC8530078	Dual pack Kit, Fits Liebel-Flarsheim Opistar Injector	96" Coiled Y-Line with 2 Check Valves and 2 Short Spikes	60 ml	50 ea/cs





DYNJHFLOW



DYNJYCLICKA



DYNJCAPMW



DYNJADAPMM



DYNJADAPMF



DYNJRADCOMP1

## Vascular Access Accessories

Designed to enhance your procedures and make your standard or custom setup compatible with your other medical components.

### Hemostatic Flow Valve

Luer-activated, self-sealing device.

» Guidewire compatibility 0.014" - 0.038" 1200 PSI

Item No.	Description	Pkg.	Nonsterile No.
DYNJHFLOW	H-Flow Valve	25 ea/cs	93956

### Hemostatic Click Device

Help prevent back bleeding in closed position and minimize back bleeding in semi-open position.

» 15 ml per second

» 600 PSI

» Guidewire compatibility 0.014 – 0.038"

Item No.	Description	Pkg.	Nonsterile No.
DYNJCLICKA	Hemostatic Click Device, Y-Click Tech	24 ea/cs	95734

### Connectors and Caps

Universal male - female fittings to ensure trusted functionality.

» Polycarbonate material

» Adaptors rated to 1200 PSI

Item No.	Description	Pkg.	Nonsterile No.
DYNJCAPMW	Male White Cap	25 ea/cs	96590
DYNJADAPMM	Adaptor M/M w/Rotator	50 ea/cs	94412
DYNJADAPMF	Adaptor M/F w/Rotator	50 ea/cs	17292

### Radial Compression Band

The solution for radial procedures that provides adequate pressure and atraumatic compression to assist with hemostasis post procedure.



Item No.	Description	Pkg.	Nonsterile No.
DYNJRADCOMP1	Radial Compression Band	100 ea/cs	122809

#### Did you know?

Medline is the market leader in procedure trays.

# Vascular Access Components for Custom or Standard Packs and Kits

Our product portfolio includes everything you need for fluid management in your Cath & Electrophysiology labs and Interventional Radiology. From manifold kits to enhanced custom packs, our team of vascular specialists works with you to tailor solutions that help improve procedural efficiency at a value.



*Products and services for all your diagnostic cardiology needs.  
For more information, contact your Medline representative.*



**Medline Industries, Inc.** Three Lakes Drive, Northfield, IL 60093

Medline United States  
1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

Medline Canada  
1-800-396-6996  
medline.ca | canada@medline.com

Medline México  
01-800-831-0898  
medlinemexico.com | mexico@medline.com

FOLLOW US

