



USP 800 protection, head to toe.

Minimize your risk to chemical exposure.

The risk is real. Are you ready for the new standard?

There's no safe level of exposure to hazardous chemotherapy drugs, but the right personal protective equipment (PPE) can minimize the risk. From face masks and gowns to gloves and footwear, count on our USP 800 recommended PPE to protect your team. So they can handle chemicals with confidence.

Read the new USP 800 guidelines at usp.org.

8 million


U.S. healthcare workers
potentially exposed to
hazardous drugs annually.¹



More than

12 billion

doses of hazardous drugs
are handled by US providers
each year.²



Outfit your staff for safety.

All recommendations below are based on USP 800.*

Gloves: Wear two pairs of disposable nitrile gloves tested for use against chemotherapy drugs per ASTM D6978.

Item Name	Item No.	Description	Sizes
FitGuard Touch Nitrile Exam Gloves	FG3002	Glove, exam, nitrile, FG Touch, 300, M	XS-XL
Restore Sense Nitrile Exam Gloves with Oatmeal	OAT2585	Glove, exam, nitrile, Restore, oat, sense, M	XS-XL
FitGuard Prime Nitrile Exam Gloves	FG2302	Glove, exam, nitrile, FitGuard Prime, M	XS-XL
Sterile Extended Cuff Nitrile Exam Gloves, Singles	MDS198515	Glove, exam, nitrile, sterile singles, M	S-XL
Sterile Extended Cuff Nitrile Exam Gloves, Pairs	MDS198415	Glove, exam, nitrile, sterile pairs, M	S-XL
Sensicare Silk Nitrile Exam Gloves	MDS2583	Glove, exam, nitrile, Sensicare Silk, M	XS-XL



For a complete list of exam gloves that are compliant with USP 800, please contact your Medline representative.

Gowns: Don a disposable, closed-back gown proven to protect against chemotherapy drugs.

Item Name	Item No.	Description	Sizes
Poly-Coated Chemotherapy Gowns	NON25457	Chemo, protective polyethylene coated, blue, L	M-XXL
Breathable Film Chemotherapy Gowns	NON27457	Chemo, premium breathable material, blue, L	M-XXL




Head and Foot Protection: Keep hair under cover and outfit closed-toe shoes with booties to reduce exposure.

Item Name	Item No.	Description	Sizes
Pro Series Bouffant Caps	NON30233B	Bouffant cap, medium-weight, blue	One size fits most
Polyethylene Shoe Covers	CRI2010	Shoe cover, polyethylene, blue	One size fits most
Premium Non-Skid Multi-Layer Shoe Covers	NON29858	Shoe cover, heavyweight SMS, blue	Regular-XL
Non-Skid Fluid Proof Boot Covers	NON27200	Boot cover, bilaminate, blue	Regular-XL



Eye and Face Protection: Wear a face shield, respirator, and goggles for splash protection.

Item Name	Item No.	Description
Fluid Protective Goggles	NON24775	Anti-fog goggles with indirect vent
¾ Length Face Shields	NONFS400	6.5" face shield
Full Length Face Shields	NONFS300	7.5" face shield
Extra-Long Face Shields	NONFS100	9.5" face shield
N95 Cone-Style Respirator Masks	NON24506A	Cone style N95 respirator, regular size
N95 Pouch Style Respirator Masks	RP88020	Pouch style N95 respirator, regular size



Wear assurance against exposure. For a free assessment, contact your Medline representative or call 1-800-MEDLINE.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US    

*Refer to section 7 of USP 800 for complete Personal Protective Equipment guidelines. Face shields alone do not provide full eye and face protection. Use in combination with goggles.

REFERENCES. 1. Centers for Disease Control: <https://www.cdc.gov/niosh/topics/hazdrug/> 2. USP 800 guidelines: <http://www.usp.org/sites/default/files/usp/document/our-work/healthcare-quality-safety/800-know-your-exposure-to-hazardous-drugs.pdf>. Both accessed May 10, 2018. © 2018 Medline Industries, Inc. All rights reserved. Medline is a registered trademark of Medline Industries, Inc. / MKT1895654 / e18189 / SG / 49