

Surgical glove ordering guide



Polyisoprene



SensiCare PI with Aloe

Item No.	MSG1055–MSG1095
Thickness	9.4 mil
Surface finish	Smooth
Color	White
Pkg.	25/bx, 100/cs
Uses	General surgery



SensiCare PI

Item No.	MSG9055–MSG9090
Thickness	9.1 mil
Surface finish	Smooth
Color	Cream
Pkg.	50/bx, 200/cs
Uses	General surgery



SensiCare PI SLT

Item No.	MSG1555–MSG1590
Thickness	9.1 mil
Surface finish	Textured
Color	Cream
Pkg.	50/bx, 200/cs
Uses	General surgery



SensiCare PI Green with Aloe

Item No.	MSG9255–MSG9290
Thickness	7.9 mil
Surface finish	Smooth
Color	Green
Pkg.	50/bx, 200/cs
Uses	Underglove



SensiCare PI Micro

Item No.	MSG9655–MSG9690
Thickness	7.9 mil
Surface finish	Smooth
Color	Cream
Pkg.	50/bx, 200/cs
Uses	Neuro, vascular, ophthalmic and cardiac procedures



SensiCare PI Ortho

Item No.	MSG9460–MSG9490
Thickness	12.6 mil
Surface finish	Smooth
Color	Cream
Pkg.	50/bx, 200/cs
Uses	Orthopedics, trauma and high-risk procedures



SensiCare PI Ortho LT*

Item No.	MSG9960–MSG9990
Thickness	12 mil
Surface finish	Textured
Color	Brown
Pkg.	50/bx, 200/cs
Uses	Orthopedics, trauma and high-risk procedures



SensiCare PI Shield**

Item No.	MSG3965–MSG3990
Thickness	11.8 mil
Surface finish	Textured
Color	Black
Pkg.	5/bx, 20/cs
Uses	Fluoroscopic procedures

Neoprene



SensiCare Neoprene

Item No.	MSG6355–MSG6390
Thickness	7.8 mil
Surface finish	Smooth
Color	Brown
Pkg.	50/bx, 200/cs
Uses	General surgery



DermAssure Green***

Item No.	MSG6555–MSG6590
Thickness	6.7 mil
Surface finish	Smooth
Color	Green
Pkg.	50/bx, 200/cs
Uses	Outer or underglove, for use with skin sensitivities

***Not made with chemical accelerators



Signature Latex Essential

Item No.	MSG5955–MSG5990
Thickness	9.1 mil
Surface finish	Smooth
Color	White
Pkg.	50/bx, 200/cs
Uses	General surgery



Triumph Micro

Item No.	MSG2355–MSG2390
Thickness	7.5 mil
Surface finish	Smooth
Color	White
Pkg.	50/bx, 200/cs
Uses	Neuro, vascular, ophthalmic and cardiac procedures

Latex

Make a
material
difference




◆ Chemo-tested per ASTM D6978-05

*Material includes Polychloroprene blend

**Material includes Tungsten blend



Make a **material** difference



Scan here
to learn more

medline.com/go/surgicalgloves

